

Product Datasheet

Human Polycomb protein SUZ12 (SUZ12) ELISA Kit (orb1656753)

Description	Human Polycomb protein SUZ12 (SUZ12) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.113 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.