

Product Datasheet

TOM1L2 antibody (orb165677)

Description	Rabbit polyclonal antibody to TOM1L2 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TOM1L2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55556
Uniprot ID	Q6ZVM7
NCBI	NP_001028723.1 , NP_001076437.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

