

Product Datasheet

Mouse Serine/threonine-protein kinase 11 (STK11) ELISA Kit (orb1656789)

Description Mouse Serine/threonine-protein kinase 11 (STK11) ELISA Kit

Reactivity Mouse

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL. Sensitivity: 0.051 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.