

## Product Datasheet

# Human Metalloreductase STEAP2 (STEAP2) ELISA Kit (orb1656792)

<b>Description</b>	Human Metalloreductase STEAP2 (STEAP2) ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.063 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids.
<b>Expiration Date</b>	Please enquire.