

Product Datasheet

Human Tomoregulin-2 (TMEFF2) ELISA Kit (orb1657022)

| | |
|----------------------------|--|
| Description | Human Tomoregulin-2 (TMEFF2) ELISA Kit |
| Reactivity | Human |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids. |
| Expiration Date | Please enquire. |