

Product Datasheet

Human Thymosin beta-10 (TMSB10) ELISA Kit (orb1657165)

Description Human Thymosin beta-10 (TMSB10)
ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: Request Information.
Sensitivity: Request Information.
Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.