

Product Datasheet

Human Tumor necrosis factor receptor superfamily member 25 (TNFRSF25) ELISA Kit (orb1657257)

Description	Human Tumor necrosis factor receptor superfamily member 25 (TNFRSF25) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.