

## Product Datasheet

# Fish Tumor Necrosis Factor $\alpha$ (TNF- $\alpha$ ) ELISA Kit (orb1657281)

<b>Description</b>	Fish Tumor Necrosis Factor $\alpha$ (TNF- $\alpha$ ) ELISA Kit
<b>Reactivity</b>	Fish
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 7.8-500 pg/mL. Sensitivity: 3.7 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids.
<b>Expiration Date</b>	Please enquire.