

Product Datasheet

Rat Transient receptor potential cation channel subfamily M member 4 (TRPM4) ELISA Kit (orb1657475)

Description	Rat Transient receptor potential cation channel subfamily M member 4 (TRPM4) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.055 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.