

## Product Datasheet

### Rat Vitamin B6 (VB6) ELISA Kit (orb1657573)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rat Vitamin B6 (VB6) ELISA Kit   |
| <b>Reactivity</b>          | Rat  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Detect Range: 3.70-300 ng/mL.<br>Sensitivity: 1.42 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.<br>Sample Volume: 1-200 $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.  |