

Product Datasheet

Human Vitamin C (VC) ELISA Kit (orb1657574)

Description	Human Vitamin C (VC) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 493.83-40000 ng/mL. Sensitivity: 192.1 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.