

Product Datasheet

Rat Vitamin D3 receptor (VDR) ELISA Kit (orb1657596)

Description Rat Vitamin D3 receptor (VDR) ELISA Kit

Reactivity Rat

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.625-40 ng/mL.
Sensitivity: 0.244 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L.
Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.