

Product Datasheet

Mouse Meprinbeta/MEP1B ELISA Kit (orb1657684)

Description	Mouse Meprinbeta/MEP1B ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156 ng/mL-10 ng/mL. Sensitivity: 10 pg/mL. Gene_name: Mep1b. Gene_full_name: meprin 1 beta. Protein_name: Meprin A subunit beta. Sample Type: cell culture supernates cell lysates tissue serum plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat.
Expiration Date	Please enquire.