

## Product Datasheet

### Mouse Meprinalpha/MEP1A ELISA Kit (orb1657685)

<b>Description</b>	Mouse Meprinalpha/MEP1A ELISA Kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.156 ng/mL-10 ng/mL. Sensitivity: 10 pg/mL. Gene_name: Mep1a. Gene_full_name: meprin 1 alpha. Protein_name: Meprin A subunit alpha. Sample Type: cell culture supernates cell lysates tissue serum plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat.
<b>Expiration Date</b>	Please enquire.