

Product Datasheet

Human GITRL/TNFSF18 ELISA Kit (orb1657689)

| | |
|----------------------------|--|
| Description | Human GITRL/TNFSF18 ELISA Kit |
| Reactivity | Human |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 125 pg/mL-8000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: TNFSF18. Gene_full_name: TNF superfamily member 18. Protein_name: Tumor necrosis factor ligand superfamily member 18. Sample Type: cell culture supernates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. |
| Expiration Date | Please enquire. |