

Product Datasheet

Human PD-L2/B7-DC ELISA Kit (orb1657704)

Description	Human PD-L2/B7-DC ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5 pg/mL-4000 pg/mL. Sensitivity: 12 pg/mL. Gene_name: PDCD1LG2. Gene_full_name: programmed cell death 1 ligand 2. Protein_name: Programmed cell death 1 ligand 2. Sample Type: cell culture supernates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat.
Expiration Date	Please enquire.