

Product Datasheet

Human TNFRSF19/TROY ELISA Kit (orb1657714)

Description	Human TNFRSF19/TROY ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5 pg/mL-4000 pg/mL. Sensitivity: 12 pg/mL. Gene_name: TNFRSF19. Gene_full_name: Tumor necrosis factor receptor superfamily member 19. Protein_name: Tumor necrosis factor receptor superfamily member 19. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Tissue_specificity: Highly expressed in prostate. Detected at lower levels in thymus spleen testis uterus small intestine colon and peripheral blood leukocytes.
Expiration Date	Please enquire.