

Product Datasheet

Mouse Legumain(total) ELISA Kit (orb1657727)

Description	Mouse Legumain(total) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: Lgmn. Gene_full_name: Legumain. Protein_name: Legumain. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Tissue_specificity: Detected in kidney proximal tubules (atprotein level). Ubiquitous. Particularly abundant in kidney andplacenta.
Expiration Date	Please enquire.