

Product Datasheet

Human Nectin-4/PVRL4 ELISA Kit (orb1657743)

Description	Human Nectin-4/PVRL4 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: NECTIN4. Gene_full_name: Nectin-4. Protein_name: Nectin-4. Sample Type: cell culture supernates serum plasma (heparin EDTA) cell lysates tissue homogenates and urine. Tmb_incubation: 15-20 min. Research_category: NECTIN4. Tissue_specificity: Predominantly expressed in placenta. Notdetected in normal breast epithelium but expressed in breastcarcinoma.
Expiration Date	Please enquire.