

Product Datasheet

Mouse R-Spondin-3 ELISA Kit (orb1657749)

Description	Mouse R-Spondin-3 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: Rspo3. Gene_full_name: R-spondin-3. Protein_name: R-spondin-3. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Research_category: stem cells signaling pathways wnt secreted developmental biology reproduction placental development. Tissue_specificity: Highly expressed in endothelial cells(PubMed:26766444).
Expiration Date	Please enquire.