

## Product Datasheet

### Human TREML1 ELISA Kit (orb1657753)

<b>Description</b>	Human TREML1 ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: TREML1. Gene_full_name: Trem-like transcript 1 protein. Protein_name: Trem-like transcript 1 protein. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Research_category: immunology innate immunity macrophage inflamm. cardiovascular blood platelets signal transduction protein phosphorylation ptyr. Tissue_specificity: Detected in platelets monocytic leukemia andin T-cell leukemia.
<b>Expiration Date</b>	Please enquire.