

## Product Datasheet

### Mouse FCRN/FCGRT ELISA Kit (orb1657789)

<b>Description</b>	Mouse FCRN/FCGRT ELISA Kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 93.7 pg/mL-6000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: Fcgrt. Gene_full_name: IgG receptor FcRn large subunit p51. Protein_name: IgG receptor FcRn large subunit p51. Sample Type: cell culture supernates serum plasma (heparin EDTA) and cell lysates. Tmb_incubation: 20-25 min. Research_category: immunology immunoglobulins receptors heavy chain igg adaptive immunity mhc class i. Tissue_specificity: Intestinal epithelium of suckling rodents.Expressed in neonatal intestine and fetal yolk sac.
<b>Expiration Date</b>	Please enquire.