

Product Datasheet

Mouse PD-L1/B7-H1 ELISA Kit (orb1657792)

Description	Mouse PD-L1/B7-H1 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: Cd274. Gene_full_name: Programmed cell death 1 ligand 1. Protein_name: Programmed cell death 1 ligand 1. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Research_category: CD274. Tissue_specificity: Highly expressed in the heart thymus skeletal muscle and lung. Weakly expressed in the kidney spleen thyroid and liver. Expressed on activated dendritic cells B-cells and macrophages. Expressed in numerous tumor cells lines of lymphoid origin.
Expiration Date	Please enquire.