

Product Datasheet

Human PD-L1/B7-H1 ELISA Kit (orb1657793)

Description	Human PD-L1/B7-H1 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5 pg/mL-4000 pg/mL. Sensitivity: 12 pg/mL. Gene_name: CD274. Gene_full_name: Programmed cell death 1 ligand 1. Protein_name: Programmed cell death 1 ligand 1. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Research_category: CD274. Tissue_specificity: Highly expressed in the heart skeletal muscle placenta and lung. Weakly expressed in the thymus spleen kidney and liver. Expressed on activated T- and B-cells dendritic cells keratinocytes and monocytes.
Expiration Date	Please enquire.