

## Product Datasheet

# Mouse Elastase/ELANE/ELA2 ELISA Kit (orb1657795)

**Description** Mouse Elastase/ELANE/ELA2 ELISA Kit

**Reactivity** Mouse

**Tested Applications** ELISA

**Note** For research use only

**Application notes** Detect Range: 31.2 pg/mL-2000 pg/mL.  
Sensitivity: 10 pg/mL. Gene\_name:  
Elane. Gene\_full\_name: Neutrophil  
elastase. Protein\_name: Neutrophil  
elastase. Sample Type: cell culture  
supernates cell lysates tissue  
homogenates serum and plasma  
(heparin EDTA). Tmb\_incubation: 20-25  
min. Research\_category: ELANE.

**Expiration Date** Please enquire.