

Product Datasheet

Rat Endostatin ELISA Kit (orb1657832)

Description	Rat Endostatin ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 15 pg/mL. Gene_name: Col18a1. Gene_full_name: collagen type XVIII alpha 1. Protein_name: Collagen alpha-1(XVIII) chain. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min.
Expiration Date	Please enquire.