

Product Datasheet

Human CRLF2/TSLP R ELISA Kit (orb1657893)

| | |
|----------------------------|--|
| Description | Human CRLF2/TSLP R ELISA Kit |
| Reactivity | Human |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: CRLF2. Gene_full_name: Cytokine receptor-like factor 2. Protein_name: Cytokine receptor-like factor 2. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity cytokines. Sequence_similarities: Belongs to the type I cytokine receptor family. Type 5 subfamily. Tissue_specificity: Expressed in heart skeletal muscle kidney and adult and fetal liver. Primarily expressed in dendrites and monocytes. Weakly expressed in T-cells. |
| Expiration Date | Please enquire. |