

## Product Datasheet

### Mouse FOLR1 ELISA Kit (orb1657915)

<b>Description</b>	Mouse FOLR1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 46.9 pg/mL-3000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: FOLR1. Gene_full_name: Folate receptor alpha. Protein_name: Folate receptor alpha. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: FOLR1. Tissue_specificity: Detected in kidney proximal tubules (atprotein level).
<b>Expiration Date</b>	Please enquire.