

Product Datasheet

Human Galectin-3BP ELISA Kit (orb1657918)

Description	Human Galectin-3BP ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5 pg/mL-4000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: LGALS3BP. Gene_full_name: Galectin-3-binding protein. Protein_name: Galectin-3-binding protein. Sample Type: cell culture supernates serum and plasma (heparin). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: cancer cell cycle cell division tumor biomarkers oncoproteins. Tissue_specificity: Ubiquitous. Detected in body fluids such as semen milk serum tears saliva and urine. Expressed by keratinocytes and fibroblasts.
Expiration Date	Please enquire.