

Product Datasheet

Mouse CCL9 ELISA Kit (orb1657928)

Description	Mouse CCL9 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	<p>Detect Range: 7.8 pg/mL-500 pg/mL. Sensitivity: 1 pg/mL. Gene_name: Ccl9. Gene_full_name: C-C motif chemokine 9. Protein_name: C-C motif chemokine 9. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity macrophage inflamm. chemokines beta chemokines (cc) cardiovascular atherosclerosis vascular inflammation leukocyte recruitment kits lysates other kits elisa kits cytokines and cytokine receptors elisa kits. Tissue_specificity: Expressed mainly in the liver lung and thymus although some expression has been detected in a widevariety of tissues except brain.</p>
Expiration Date	Please enquire.