

Product Datasheet

Rat MIP-1Alpha/CCL3 ELISA Kit (orb1657933)

Description	Rat MIP-1Alpha/CCL3 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.8 pg/mL-500 pg/mL. Sensitivity: 1 pg/mL. Gene_name: CCL3. Gene_full_name: C-C motif chemokine 3. Protein_name: C-C motif chemokine 3. Sample Type: cell culture supernates serum and plasma (heparin EDTA citrate). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity macrophage inflamm. chemokines beta chemokines (cc) kits lysates other kits elisa kits blood cell antigens elisa kits immune system diseases antiviral signaling hiv-related.
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)