

## Product Datasheet

### Mouse Progranulin ELISA Kit (orb1657948)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Mouse Progranulin ELISA Kit   |
| <b>Reactivity</b>          | Mouse   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Detect Range: 0.78 ng/mL-50 ng/mL.<br>Sensitivity: 10 pg/mL. Gene_name: GRN.<br>Gene_full_name: Granulins. Protein_name: Granulins. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: GRN. Tissue_specificity: Ubiquitous. |
| <b>Expiration Date</b>     | Please enquire.   |