

Product Datasheet

Mouse IL-17RC ELISA Kit (orb1657957)

Description	Mouse IL-17RC ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: IL17RC. Gene_full_name: Interleukin-17 receptor C. Protein_name: Interleukin-17 receptor C. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity cytokines interleukins.
Expiration Date	Please enquire.