

Product Datasheet

Human KLK12 ELISA Kit (orb1657967)

Description	Human KLK12 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: KLK12. Gene_full_name: Kallikrein-12. Protein_name: Kallikrein-12. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: signal transduction cytoskeleton ecm extracellular matrix ecm enzymes kallikreins cancer invasion microenvironment cell biology proteolysis ubiquitin proteolytic enzymes serine protease.
Expiration Date	Please enquire.