

Product Datasheet

Mouse CCL8/MCP-2 ELISA Kit (orb1657979)

| | |
|----------------------------|---|
| Description | Mouse CCL8/MCP-2 ELISA Kit |
| Reactivity | Mouse |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: CCL8. Gene_full_name: C-C motif chemokine 8. Protein_name: C-C motif chemokine 8. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity macrophage inflamm. chemokines beta chemokines (cc) kits lysates other kits elisa kits cytokines and cytokine receptors elisa kits. |
| Expiration Date | Please enquire. |