

Product Datasheet

Mouse CCL6/C10 ELISA Kit (orb1657985)

Description	Mouse CCL6/C10 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: Ccl6. Gene_full_name: C-C motif chemokine 6. Protein_name: C-C motif chemokine 6. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity chemokines beta chemokines (cc). Tissue_specificity: Expressed in myelopoietic bone marrow culturesstimulated by GM-CSF.
Expiration Date	Please enquire.