

## Product Datasheet

### Human Progranulin ELISA Kit (orb1658021)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human Progranulin ELISA Kit   |
| <b>Reactivity</b>          | Human   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Detect Range: 1.56 ng/mL-100 ng/mL.<br>Sensitivity: 50 pg/mL. Gene_name: GRN.<br>Gene_full_name: Granulins. Protein_name: Granulins. Sample Type: cell culture supernates serum plasma (heparin EDTA) saliva urine and human milk.<br>Tmb_incubation: 25-30 min.<br>Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat.<br>Research_category: GRN.<br>Sequence_similarities: Belongs to the granulin family. Tissue_specificity: In myelogenous leukemic cell lines of promonocytic promyelocytic and proerythroid lineage in fibroblasts and very strongly in epithelial cell lines. Present in inflammatory cells and bone marrow. Highest levels in kidney. |
| <b>Expiration Date</b>     | Please enquire.   |