

Product Datasheet

Rat CCL4/MIP-1 Beta ELISA Kit (orb1658046)

Description	Rat CCL4/MIP-1 Beta ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: CCL4. Gene_full_name: C-C motif chemokine 4. Protein_name: C-C motif chemokine 4. Sample Type: cell culture supernates serum and plasma (heparin EDTA citrate). Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity macrophage inflamm. chemokines beta chemokines (cc) kits lysates other kits elisa kits cytokines and cytokine receptors elisa kits immune system diseases antiviral signaling hiv-related.
Expiration Date	Please enquire.