

Product Datasheet

Human Angiostatin K1-3 ELISA Kit (orb1658064)

Description	Human Angiostatin K1-3 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: PLG. Gene_full_name: Plasminogen. Protein_name: Plasminogen. Sample Type: serum plasma and cell culture supernates. Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: cardiovascular blood fibrinolysis thrombolysis platelets coagulation regulatory cell biology proteolysis ubiquitin proteolytic enzymes proteases protease inhibitors protease inhibitors kits lysates other kits elisa kits cancer proteins elisa kits protease inhibitors elisa kits. Sequence_similarities: Belongs to the peptidase S1 family. Plasminogen subfamily. Tissue_specificity: Present in plasma and many other extracellular fluids. It is synthesized in the liver.
Expiration Date	Please enquire.