

## Product Datasheet

### Human NOV/CCN3 ELISA Kit (orb1658096)

<b>Description</b>	Human NOV/CCN3 ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 2 pg/mL. Gene_name: NOV. Gene_full_name: Protein NOV homolog. Protein_name: Protein NOV homolog. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: Nothing Found. Sequence_similarities: Belongs to the CCN family. Tissue_specificity: Expressed in bone marrow thymic cells and nephroblastoma. Increased expression in Wilms tumor of the stromal type.
<b>Expiration Date</b>	Please enquire.