

Product Datasheet

Mouse LOX-1/OLR1 ELISA Kit (orb1658103)

| | |
|----------------------------|--|
| Description | Mouse LOX-1/OLR1 ELISA Kit |
| Reactivity | Mouse |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 1 pg/mL. Gene_name: OLR1. Gene_full_name: Oxidized low-density lipoprotein receptor 1. Protein_name: Oxidized low-density lipoprotein receptor 1. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: cardiovascular lipids lipoproteins fatty acids binding proteins signal transduction metabolism lipid metabolism cell biology antibodies oxidative stress kits lysates other kits elisa kits lipoprotein elisa kits. Sequence_similarities: Contains 1 C-type lectin domain. |
| Expiration Date | Please enquire. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)