

Product Datasheet

Mouse LOX-1/OLR1 ELISA Kit (orb1658103)

Description	Mouse LOX-1/OLR1 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 1 pg/mL. Gene_name: OLR1. Gene_full_name: Oxidized low-density lipoprotein receptor 1. Protein_name: Oxidized low-density lipoprotein receptor 1. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: cardiovascular lipids lipoproteins fatty acids binding proteins signal transduction metabolism lipid metabolism cell biology antibodies oxidative stress kits lysates other kits elisa kits lipoprotein elisa kits. Sequence_similarities: Contains 1 C-type lectin domain.
Expiration Date	Please enquire.