

Product Datasheet

Rat CXCL5/LIX/ENA-78 ELISA Kit (orb1658105)

Description	Rat CXCL5/LIX/ENA-78 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 1 pg/mL. Gene_name: Cxcl5. Gene_full_name: C-X-C motif chemokine 5. Protein_name: C-X-C motif chemokine 5. Sample Type: cell culture supernates serum plasma (heparin EDTA citrate) and tissue homogenates. Tmb_incubation: 15- 20 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: Nothing Found.
Expiration Date	Please enquire.