

Product Datasheet

Mouse Mer ELISA Kit (orb1658111)

Description	Mouse Mer ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78 pg/mL-5000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: MERTK. Gene_full_name: Tyrosine-protein kinase Mer. Protein_name: Tyrosine-protein kinase Mer. Sample Type: cell culture supernates cell lysates tissue homogenates and serum. Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: neuroscience neurotransmission intracellular signaling kinases signal transduction protein phosphorylation tyrosine kinases receptor tyrosine kinases sensory system visual system. Sequence_similarities: Belongs to the protein kinase superfamily. Tyr protein kinase family. AXL UFO subfamily. Tissue_specificity: Expressed predominantly in the hematopoieticlineages: macrophages NK cells NKT cells dendritic cells andplatelets.
Expiration Date	Please enquire.