

Product Datasheet

Human IL-17RA ELISA Kit (orb1658127)

Description	Human IL-17RA ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: IL17RA. Gene_full_name: Interleukin-17 receptor A. Protein_name: Interleukin-17 receptor A. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: IL17RA. Tissue_specificity: Widely expressed.
Expiration Date	Please enquire.