

Product Datasheet

Mouse CD40/TNFRSF5 ELISA Kit (orb1658181)

Description	Mouse CD40/TNFRSF5 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 1 pg/mL. Gene_name: CD40. Gene_full_name: Tumor necrosis factor receptor superfamily member 5. Protein_name: Tumor necrosis factor receptor superfamily member 5. Sample Type: cell culture supernates cell lysates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology adaptive immunity b cells cd cardiovascular atherosclerosis vascular inflammation innate and adaptive immunity stem cells hematopoietic progenitors lymphoid b lymphocytic lineage myeloid dendritic cell lineage.
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)