

Product Datasheet

Mouse CCL28 ELISA Kit (orb1658190)

Description	Mouse CCL28 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: CCL28. Gene_full_name: C-C motif chemokine 28. Protein_name: C-C motif chemokine 28. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity chemokines beta chemokines (cc) cytokines metabolism types of disease obesity. Tissue_specificity: Mainly expressed in testis epithelial cells of normal colon kidney Peyer patches lymph nodes. Also found in lower levels in brain spleen and lung.
Expiration Date	Please enquire.