

## Product Datasheet

### Mouse TNFRSF17/BCMA ELISA Kit (orb1658202)

<b>Description</b>	Mouse TNFRSF17/BCMA ELISA Kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: TNFRSF17. Gene_full_name: Tumor necrosis factor receptor superfamily member 17. Protein_name: Tumor necrosis factor receptor superfamily member 17. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: TNFRSF17. Tissue_specificity: Detected in spleen thymus bone marrow and heart and at lower levels in kidney and lung.
<b>Expiration Date</b>	Please enquire.