

Product Datasheet

Human TNFRSF17/BCMA ELISA Kit (orb1658203)

Description	Human TNFRSF17/BCMA ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: TNFRSF17. Gene_full_name: Tumor necrosis factor receptor superfamily member 17. Protein_name: Tumor necrosis factor receptor superfamily member 17. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: TNFRSF17. Sequence_similarities: Contains 1 TNFR-Cys repeat. Tissue_specificity: Expressed in mature B-cells but not in T-cells or monocytes.
Expiration Date	Please enquire.