

Product Datasheet

Mouse Gp130/IL6ST ELISA Kit (orb1658220)

Description	Mouse Gp130/IL6ST ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 125 pg/mL-8000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: IL6ST. Gene_full_name: Interleukin-6 receptor subunit beta. Protein_name: Interleukin-6 receptor subunit beta. Sample Type: cell culture supernates serum and plasma (heparin EDTA citrate). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: cell biology apoptosis intracellular associated proteins immunology cell type markers cd cytokines signal transduction adapters transmembrane cancer cardiovascular atherosclerosis heart hypertrophy. Sequence_similarities: Belongs to the type I cytokine receptor family. Type 2 subfamily. Tissue_specificity: Found in tissues such as brain heart thymus spleen kidney lung and liver. Found in all the cell lines tested except BaF-B03. Expression not restricted to IL6-responsive cells.
Expiration Date	Please enquire.